# TympoVent<sup>®</sup>

## **Disclaimer**

Atos Medical offers no warranty - neither expressed nor implied - to the purchaser hereunder as to the lifetime of the product delivered, which may vary with individual use and biological conditions. Furthermore, Atos Medical offers no warranty of merchantability or fitness of the product for any particular purpose.

## **Patents and trademarks**

TympoVent® is a registered trademark owned by Atos Medical AB, Sweden. For information about protective rights (e.g. patents), please refer to our web page [www.atosmedical.com/patents](http://www.atosmedical.com/patents).

Manufacturer; Hersteller; Fabrikant; Fabricant; Produttore;  
 Fabricante; Fabricante; Tillverkare; Producent; Producent;  
 Valmistaja; Framleiðandi; Tootja; Ražotājs; Gamintojas;  
 Výrobce; Gyártó; Výrobca; Proizvajalec; Producent;  
 Producător; Proizvođač; Proizvođač; Κατασκευαστής;  
 Производител; Üretici; მწარმოებელი; Ἀρτοποιητής;  
 İstehsalçı; Производитель; Produzen; Pengilang; 製造業  
 者; 제조사; 製造商; 制造商;

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## Ventilation (Tympanostomy) Tubes

### Device Description:

TympoVent ventilation (tympanostomy) tubes are small tubular implants, which are placed through an incision made in the tympanic membrane in order to connect the middle ear cavity with the external auditory canal. They are used for temporary ventilation of the middle ear and/or for drainage of fluids from the middle ear into the external auditory canal.

The insertion of the device is to be conducted by a suitably qualified clinician in a hospital.

### Indications:

Pathologic conditions where ventilation tubes are indicated include for example:

- Secretory otitis media with effusion
- Recurrent acute otitis media not responding satisfactorily to conventional therapy.

### Contraindications:

- Ongoing middle ear infection
- When adverse effects may be expected (see below).

### Warnings:

Since the insertion of ventilation tubes creates an open passageway between the middle ear and the external environment, the middle ear may develop subsequent secondary infections from either water or airborne pathogens. Care should be taken to avoid fluid entering the ear to reduce the risk of exposure to waterborne pathogens.

Re-use and re-processing may cause cross-contamination and damage to the device, which could cause patient harm.

### Possible adverse effects and complications:

- Clogging of the tube; may require removal and replacement.
- Permanent perforation after the tube has extruded or been removed.
- Discharge caused by secondary infection of pathogens entering through the tube.
- Tissue reaction caused by allergy or sensitivity to certain materials.
- Damage to the ossicular chain, facial nerve, chorda tympani and a possible high jugular bulb caused by incorrect placement of the tube.
- Myringosclerotic changes shown as white plaques.
- Cholesteatoma formation caused by ingrowth of squamous epithelium along the tube into the middle ear.

## Disposal

Always follow medical practice and national requirements regarding biohazard when disposing of a used medical device.

## Device lifetime

Tympanostomy tubes (TT) are commonly in-situ between a few days (early extrusion) up to 2 years and more. The tubes usually extrude spontaneously due to growth of the tympanic membrane and only in rare cases is a surgical removal needed, for example in case of ongoing otitis media. Besides patient related variables, the time the tubes are in situ also depends on:

- Type of tube: short (e.g. Donaldson, Shephard) or long (e.g. Armstrong, Straight Tube)
- Type of tube material: silicone, fluoroplastic, titanium, gold/silver.

## Reporting

Please note that any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the national authority of the country in which the user and/or patient resides.

**Atos Medical does not recommend re-sterilisation of its ventilation tubes.**

**Atos Medical assumes no liability for products, which have been re-sterilised.**

### Caution!

**TympoVent tubes containing gold and titanium are not considered safe for magnetic resonance imaging (MRI), and should be removed prior to procedure. Tubes containing only silicone or plastic are safe to use during MRI, and can remain in situ.**

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