PROVOX[®]


Coming Home[®]

Provox Luna Adhesive

Figure 1


1.1


1.2


1.3


1.4


1.5


1.6


1.7


1.8

1.9

1.10

1.11




Figure 2


2.1


2.2


2.3


2.4


2.5


2.6

2.7

2.8

2.9



Figure 3


3.1


3.2

3.3

3.4

3.5

3.6





Figure 4
4.1


4.2


4.3


4.4

4.5

4.6

4.7




Figure 4


4.8

4.9

4.10

Disclaimer

Atos Medical offers no warranty - neither expressed nor implied - to the purchaser hereunder as to the lifetime of the product delivered, which may vary with individual use and biological conditions. Furthermore, Atos Medical offers no warranty of merchantability or fitness of the product for any particular purpose.

Patents and trademarks

Provox[®], Provox[®] Luna[®] and Coming Home[®] are registered trademarks owned by Atos Medical AB, Sweden.

Micron HME[™], FlexiDerm[™], OptiDerm[™], StabiliBase[™], XtraHME[™], XtraMoist[™], XtraFlow[™] and ShowerAid[™] are trademarks of Atos Medical AB.

For information about protective rights (e.g. patents), please refer to our web page www.atosmedical.com/patents.

Manufacturer; Hersteller; Fabrikant; Fabricant; Produttore; Fabricante; Fabricante; Tillverkare; Producent; Producent; Valmistaja; Framleiðandi; Tootja; Ražotājs; Gamintojas; Výrobce; Gyártó; Výrobca; Proizvajalec; Producent; Producător; Proizvođač; Proizvođač; Κατασκευαστής; Производител; Üretici; მწარმოებელი; Ἀρτοποιητής; İstehsalçı; Производитель; Produ-
dusen; Pengilang; 製造業者; 제조사; 製造商; 制造商; היהצור; الجهة المصنعة

Date of manufacture; Herstellungsdatum; Productiedatum; Date de fabrication; Data di fabbricazione; Fecha de fabricación; Data de fabrico; Tillverkningsdatum; Fremstillingsdato; Produksjonsdato ;Valmistuspäivä; Framleiðsludagsetning; Valmistamise kuupäev; Ražošanas datums; Pagaminimo data; Datum výroby; Gyártás időpontja; Dátum výroby; Datum proizvodnje; Data produkcji; Data fabricației; Datum proizvodnje; Datum proizvodnje; Ημερομηνία κατασκευής; Дата на производство; Üretim tarihi; წარმოების თარიღი; Ἀρτοποιημάτων ἀνάσπαση; İstehsal tarixi; Дата изготовления; Tanggal produksi; Tarikh pembuatan; 製造日; 제조 일자; 製造日期; 製造日期; תאריך ייצור; تاريخ التصنيع

Use-by date; Verfallsdatum; Houdbaarheidsdatum; Date limite d'utilisation; Data di scadenza; Fecha de caducidad; Data de validade; Används senast; Anvendes inden; Brukes innen; Viimeinen käyttöpäivä; Sīdasti notkunardagur; Kõlblik kuni; Derīguma termiņš; Naudoti iki; Datum expirace; Felhasználatóság dátuma;

Do not re-use; Nur zum Einmalgebrauch; Niet hergebruiken; Ne pas réutiliser; Non riutilizzare; No reutilizar; Não reutilizar; Återanvänd ej; Må ikke genanvendes; Skal ikke gjenbrukes; Ei uudelleenkäytettävä; Má ekki endurnýta; Mitte taaskasutada; Neizmantot atkārtoti; Pakartotinai nenaudoti; Nepoužívejte opakovaně; Ne használja fel újra; Nepoužívajte opakovane; Ne uporabljajte ponovno; Nie używać ponownie; A nu se reutiliza; Nemojte ponovno upotrebljavati; Ne koristiti ponovo; Μην επαναχρησιμοποιείτε; Не използвайте повторно; Tekrar kullanmayın; არ გამოიყენოთ ხელახლა; Մեկանգամյա օգտագործման համար; Yenidən istifadə etməyin; Не подлежит повторному использованию; Jangan gunakan ulang; Jangan guna semula; 再利用不可; 재사용하지 마십시오; 請勿重複使用; 请勿重复使用;

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Intended use

Provox® Coming Home® is an assortment of products and information for newly laryngectomized patients. It provides guidance on product use, lungrehabilitation, and stoma care at home.

How to use Provox® Coming Home®

Provox Coming Home is intended for single patient use.

The products included are selected to be used every day after a total laryngectomy. See instructions on how to use the different products below. Additional information can be found in Provox Coming Home Booklet.

Disposal

Always follow medical practice and national requirements regarding biohazard when disposing of a used medical device.

User assistance information

For additional help or information, please see Contact card provided.

Different products for different situations

In Provox Coming home you will find products designed to help in different situations; Shower, Care, Day, Night, Going Out and Moving On. The products can be found in the different pockets listed below.

Reporting

Please note that any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the national authority of the country in which the user and/or patient resides.

Table 1

Product	Pocket					
	Shower	Care	Day	Night	Going Out	Moving On
Provox ShowerAid	X					
Provox Luna ShowerAid	X					
Provox Cleaning Towel		X				
Provox Skin Barrier		X				
Provox Adhesive Remover		X				
Provox FlexiDerm			X		X	
Provox OptiDerm			X			
Provox XtraFlow			X			
Provox Luna				X		
Provox Micron					X	
Provox StabiliBase						X

Provox® ShowerAid™

Intended use

Provox ShowerAid is used to temporarily replace Provox XtraMoist/Provox XtraFlow/Provox FreeHands/Provox Micron HME during showering. The device is for single patient use.

Provox® Luna® ShowerAid

Intended use

Provox Luna ShowerAid is used with Provox Luna Adhesive while taking a shower to avoid water from entering the stoma. Single patient use.

How to use

Attach and remove Provox ShowerAid

Prior to entering the shower, remove the HME and insert Provox ShowerAid with the opening facing down (Fig. 1.1-1.2). After showering, remove Provox ShowerAid and insert an HME (Fig. 1.3 – 1.4).

Attach and remove Provox Luna ShowerAid

Prior to entering the shower, remove Provox Luna HME, attach Provox Adhesive strip (Fig. 1.5-1.6) and insert Provox Luna ShowerAid with the opening facing down (Fig. 1.1-1.2). After showering remove Provox Luna ShowerAid (Fig 1.3).

Cleaning and disinfection

Clean the device after each use (Fig. 1.7-1.9). Disinfect monthly using 70% Ethanol or 70 % Isopropylalcohol for 10 minutes or 3% Hydrogenperoxide for 60 minutes (Fig. 1.10 – 1.11). When the device shows signs of damage, it must be discarded. Replace Provox ShowerAid or Provox Luna ShowerAid at least yearly.

WARNING

- Single patient use only. Reuse between patients may cause cross-contamination.

Provox® Cleaning Towel

Intended use

Provox Cleaning Towel is intended for cleaning around the stoma, it will remove oil from the skin. They are intended to use before application of Provox Adhesives.

How to use Provox Cleaning Towel

Application

Clean the skin around the stoma with Provox Cleaning Towel and let dry.

Provox® Skin Barrier

Intended use

Provox Skin Barrier is a single use wipe for laryngectomized patients that forms a barrier between Provox Adhesive and the skin.

How to use Provox Skin Barrier

Application:

Clean the skin around the stoma and let dry. Apply Provox Skin Barrier to the skin around the stoma that will be covered by the adhesive. Allow the barrier to dry. Apply a Provox Adhesive.

CONTRAINDICATION

Hypersensitive or breached skin or known allergies to ingredients (Hexamethyldisiloxane, Isopropyl myristate, Trimethylsiloxysilicate)

WARNINGS

- For external use only.
- Keep out of reach of children.
- Avoid contact with eyes. In the case of accidental contact, flush eyes well with water.
- Do not apply directly to open wounds and mucous membranes.
- Single patient use only. Re-use may cause cross contamination.
- Do not cut the tissue in pieces, small pieces may be aspirated

PRECAUTIONS

- Avoid deep inhalation of solvent during application and be careful that liquid does not drip into stoma.
- In case you experience signs of irritation or redness, discontinue use of product. If the irritation persists, consult your physician.

Provox® Adhesive Remover

Intended use

Provox Adhesive Remover is a single use wipe to help laryngectomized patients remove Provox Adhesives and Provox Silicone Glue.

How to use Provox Adhesive Remover

Application

When using a Provox FlexiDerm or Provox StabiliBase, first apply Provox Adhesive Remover on top of the adhesive. This is not necessary for the other adhesives. Then, for all Provox adhesives, grasp the finger lift tab and apply Provox Adhesive Remover on the edge and underneath the adhesive. Carefully remove the adhesive and adhesive residue from the skin. Clean the skin with water and soap or use Provox Cleaning Towel afterwards.

CONTRAINDICATION

Hypersensitive or breached skin or known allergies to ingredients (Hexamethyldisiloxane, Isopropyl myristate).

WARNINGS

- For external use only.
- Keep out of reach of children.
- Avoid contact with eyes. In the case of accidental contact, flush eyes well with water.
- Do not apply directly to open wounds and mucous membranes.
- Single patient use only. Re-use may cause cross contamination.
- Do not cut the tissue in pieces, small pieces may be aspirated.

PRECAUTIONS

- Avoid deep inhalation of solvent during application and be careful that liquid does not drip into stoma.
- In case you experience signs of irritation or redness, discontinue use of product. If the irritation persists, consult your physician.

Provox® Adhesives

Provox® FlexiDerm™, Provox® OptiDerm™, Provox® StabiliBase™

Intended use

Provox Adhesives are single use devices intended for laryngectomized patients breathing through a tracheostoma. The devices are attached to the skin around the tracheostoma in order to provide attachment of components of Provox HME System.

How to attach Provox FlexiDerm

To apply the adhesive to your skin, first make sure the skin is clean and dry. If needed, you can apply skin protection products (described above). Then remove the back liner and apply Provox FlexiDerm around your stoma. Gently massage the adhesive onto the skin to improve adherence (Fig. 2.2 – 2.4).

How to attach Provox OptiDerm

To apply the adhesive to your skin, first make sure the skin is clean and dry. Pre-warm the adhesive between your hands to improve the adhesion. Then remove the back liner and apply Provox OptiDerm around your stoma. Gently massage the adhesive onto the skin to improve adherence (Fig. 2.1 – 2.4). Remove Provox OptiDerm when it is loose or dirty. Carefully remove the adhesive from the skin, using the finger lift tab.

Caution: When using Provox OptiDerm in the post-operative period or on sensitive skin, it should be removed very slowly and carefully. When cleaning the skin from e.g. residual glue, prevent particles/fluids from entering your stoma.

How to attach Provox StabiliBase

To apply the adhesive to your skin, first make sure the skin is clean and dry. If needed, you can apply skin protection products (described above). Then remove the center piece of the back-liner (Fig. 2.5-2.6). Apply the adhesive to the skin around your stoma (Fig. 2.7), and then remove the side pieces of the back liner (Fig. 2.8). Gently massage the adhesive onto the skin to improve adherence (Fig. 2.9).

How to remove Provox Adhesives

Remove Provox Adhesive when it has become loose or dirty. Carefully remove the adhesive from the skin, using the finger lift tab. An adhesive remover (described above) may be helpful for removing adhesives or glue. Always clean the skin with a Provox Cleaning Towel and/or soap and water after use of an adhesive remover. Dry the area carefully.

Caution: When cleaning the skin from e.g. residual glue, prevent particles/fluids from entering your stoma.

WARNINGS

- The adhesive may irritate the skin. Stop using the adhesive if skin irritation develops and consult your clinician.
- Do not use Provox Adhesives during radiotherapy that is delivered to the area covered by the adhesive. Consult your clinician before you resume use of the adhesive after radiation therapy.
- Reuse of the adhesive, by yourself or someone else, may cause transfer of micro-organisms which can lead to infections.
- Reuse may also cause reduced efficiency of the adhesive properties which can lead to air leakage during speech and reduced efficiency of the attached HME.
- Only use genuine Provox system components that are intended for use with Provox Adhesives. Other devices may cause personal injury or damage to the products.

Provox® XtraFlow™

Intended use

Provox XtraHME Cassette is a single use, specialized device intended for patients breathing through a tracheostoma. It is a heat and moisture exchanger (HME) that heats and humidifies inhaled air by retaining heat and moisture from exhaled air in the device. It partially restores lost breathing resistance. For patients with a voice prosthesis or surgical speech fistula it may also facilitate voicing.

How to use Provox XtraFlow

Attach and remove

Place the HME in the holder of your attachment device (e.g. Provox Adhesive or Provox LaryTube) and breathe normally (Fig. 3.1).

To remove the HME, hold the attachment device in place with two fingers and remove the HME from the holder (Fig. 3.3).

Speaking with a voice prosthesis

If you have a voice prosthesis and you have been cleared to use it for speaking, the lid of Provox XtraFlow can be pressed down with your finger to occlude the stoma for speaking. After releasing the finger pressure, the lid opens and you can breathe (Fig. 3.2).

CONTRAINDICATIONS

- This device shall not be used by patients who are unable to handle or remove the device themselves when needed, unless the patient is under constant supervision of a clinician or a trained caregiver. For example: patients who are unable to move their arms, patients with decreased levels of consciousness, or patients with diseases that put them at a risk for unpredictable periodic loss of consciousness.
- Do not use on patients with a low tidal volume, as the added dead space (5 ml) may cause CO₂ (carbon dioxide) retention at too low tidal volume.

WARNING

- Be careful not to exert pressure on the lid of the HME unintentionally. Unintentional closing of the lid may cause difficulty in breathing.
- Always inform the patient, caretakers and others about the closing feature of the HME cassette to ensure that they understand its function. Closing the airway to allow voicing is a well-known feature for the laryngectomized patient with a voice prosthesis. For patients without a voice prosthesis or tracheostomized patients this feature might be unknown.

PRECAUTIONS

- Always test the function of Provox XtraFlow prior to use. The top lid should immediately return to its open position after releasing the pressure.

- Do not disassemble Provox XtraFlow since this will interfere with its proper function.
- Do not reuse Provox XtraFlow or attempt to rinse it with water or any other substance. This will substantially reduce the function of the HME. It also increases the risk of potential infections since e.g. bacteria may start to grow in the foam.
- Do not use Provox XtraFlow longer than 24 hours. This will increase the risk of potential infections since e.g. bacteria may start to grow in the foam.
- Do not administer medicated nebulizer treatment over the device since the medication can become deposited in the device.
- Do not use humidifiers or heated humidified oxygen via a mask over the tracheostoma while using the device. The HME will become too wet. If oxygen therapy is required, use only non-heated humidified oxygen.

Provox® Micron HME™

Provox Micron HME is not available in all countries. In the US this product is only available on prescription and in Japan it is not available at all.

Intended use

Provox Micron HME is a heat and moisture exchanger (HME) and air filtration device for patients breathing through a tracheostoma. Provox Micron HME partially restores lost breathing resistance. For patients with a voice prosthesis or surgical speech fistula it may also facilitate voicing.

Provox Micron HME is intended to be used with the attachment devices in Provox HME System.

How to use Provox Micron HME

Attach and remove

Provox Micron HME can easily be attached (Fig. 3.4) and removed (Fig. 3.6) from Provox HME System attachment devices when needed.

Provox Micron HME helps to filter inhaled and exhaled air through consistent normal use. Thereby, small airborne particles, e.g., bacteria, viruses, dust and pollen are restricted from passing through the device into the lungs upon inspiration (see technical data below), while also protecting others nearby during expiration.

Note: Provox Micron is not intended to be used as a Personal Protective Equipment during work that requires breathing protection.

Speaking with a voice prosthesis

If you have a voice prosthesis and you have been cleared to use it for speaking, the lid of Provox Micron HME can be pressed down with your finger to occlude the stoma for speaking. After releasing the finger pressure, the lid opens and you can breathe (Fig. 3.5).

CONTRAINDICATIONS

This device shall not be used by patients who are unable to handle or remove the device themselves when needed, unless the patient is under constant supervision of a clinician or a trained caregiver. For example: patients who are unable to move their arms, patients with decreased levels of consciousness, or patients with diseases that put them at a risk for unpredictable periodic loss of consciousness.

WARNING

Be careful not to exert pressure on the lid of the HME unintentionally. Unintentional closing of the lid may cause difficulty in breathing.

PRECAUTIONS

- Provox Micron HME provides good protection through consistent normal use, as long as there is no air leakage. However, since there are other pathways for e.g. viruses and bacteria to enter the human body, total protection can never be guaranteed.
- The same device must not be used for more than 24 hours after initial use. This can increase the risk for infection due to growth of e.g. bacteria.
- Do not wash and re-use the device. Washing the HME impairs the filtering and HME functions.
- Do not disassemble Provox Micron HME. Disassembly will destroy its function.
- Do not administer medicated nebulizer treatment over the device since the medication can be deposited in the device.
- Do not use humidifiers or heated humidified oxygen over the device since the HME will become too wet.
- Replace Provox Micron HME when needed. To ensure its proper function, the same device must not be used for more than 24 hours after initial use.

Provox® Luna®

Intended use

Provox Luna HME is a single use heat and moisture exchanger, attachable to Provox Luna Adhesive, for night-time use after total laryngectomy.

Provox Luna Adhesive is a skin-friendly, single use adhesive that provides attachment for Provox Luna HME for night-time use after total laryngectomy.

How to use Provox Luna HME and Provox Luna Adhesive

Attach and remove

Clean the skin around the stoma and let it dry before attaching Provox Luna Adhesive.

For instructions on how to apply Provox Luna adhesive, see fig 4.1-4.7.

Caution: Do not use skin barrier products or skin cream.

Place Provox Luna HME in Provox Luna Adhesive according to figure 4.8. In order to get a seal for speech, occlude Provox Luna HME by covering the openings on the sides (fig. 4.9).

Provox Luna HME can be removed (as in fig. 4.10) if mucus needs to be cleaned from the stoma. In the morning, carefully remove Provox Luna Adhesive. Using a little bit of water at the side and underneath the adhesive helps to remove the adhesive.

Caution: When taking a shower with Provox Luna Adhesive use Provox Luna ShowerAid and Provox Adhesive Strip. The hydrogel adhesive material may absorb water and detach.

CONTRAINDICATION

- The product shall not be used by patients with a decreased level of consciousness, patients with reduced mobility of the arms and/or hands, or patients who are unable to remove the device themselves.
- The product shall not be used by patients with a low tidal volume, as the added dead space (7 ml) may cause CO₂ (Carbon dioxide) retention.

WARNING

- Do not use Provox Luna HME or Provox Luna ShowerAid with any other adhesives or other accessories unless indicated in this Instructions for Use.
- Do not use Provox Luna Adhesive with any other HME or accessories unless indicated in this Instructions for Use.

- Provox Luna Adhesive should not be worn during radiotherapy treatment sessions as it will alter the delivery of radiotherapy to that area. As Provox Luna Adhesive is a hydrogel, it is possible to use this adhesive during the night between radiotherapy sessions; however, clinicians should be aware that this may not be appropriate for all patients and use of Provox Luna Adhesives needs to be assessed both on an individual basis and reviewed regularly through the course of radiotherapy treatment.

PRECAUTION

- Be careful not to exert pressure on Provox Luna HME unintentionally. Pressure may reduce the size of the side openings and cause difficulty breathing.
- Do not reuse Provox Luna HME as the risk of potential infections may increase due to bacterial colonization of the HME.
- Do not rinse Provox Luna HME with water or any other substance since this will substantially reduce the function of the HME.
- Do not remove the HME foam, as the HME function will be lost.
- Provox Luna Adhesive may be used on irritated skin. However, if skin irritation persists for more than 14 days, stop using the adhesive and consult your clinician.
- If new skin irritation develops, stop using Provox Luna Adhesive and consult with your clinician.

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*Caution: United States Federal law restricts the sale, distribution and use by, or by order of a physician or a licensed practitioner of devices identified herein with the symbol "Rx". The prescription-free availability of these products outside the United States may vary from country to country.

THIS PAGE IS DELIBERATELY LEFT EMPTY

THIS PAGE IS DELIBERATELY LEFT EMPTY

THIS PAGE IS DELIBERATELY LEFT EMPTY

No. 11299, 2020-10-07

Atos
atosmedical.com

Atos Medical AB
Kraftgatan 8, SE-242 35 Hörby, Sweden
Tel: +46 (0)415 198 00 • info@atosmedical.com